

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium tetracyanopalladate(II) hydrate - 98%

Product Number: **334553** $K_2Pd(CN)_4 \cdot xH_2O$

MDL: MFCD00150432
Formula: $C_4K_2N_4Pd \cdot xH_2O$
Formula Weight: 288.69 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Crystals
ICP Major Analysis Confirms K and Pd Components	Confirmed
Carbon	13.0 - 16.1 %
Nitrogen	15.1 - 18.7 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.