

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum on alumina - extent of labeling: 1 wt. % loading, 3.2 mm, pellets

Product Number:

232114

MDL:

MFCD00011179

Pt**TEST****Specification**

Appearance (Color)

Grey to Dark Grey

Appearance (Form)

Pellets

Platinum (Pt)

0.5 - 1.5 %

Specification: PRD.1.ZQ5.10000015316

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.