


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Amino-2,3-dihydrobenzoic acid hydrochloride - ≥98%

Product Number: A3539
CAS Number: 59556-17-1
MDL: MFCD00036983
Formula: C₇H₉NO₂ · HCl
Formula Weight: 175.61 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
50 mg/ml , H ₂ O	
Solubility (Turbidity)	Clear
Water (by Karl Fischer)	≤ 10 %
Purity (Tit. by NaOH) anhyd.	≥ 98 %
Purity (TLC)	≥ 98 %
Recommended Retest Period	-----
4 years	

Specification: PRD.0.ZQ5.10000074350

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.