

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Tumor Necrosis Factor- α antibody produced in mouse - clone 28401, purified immunoglobulin, lyophilized powderProduct Number: **T6817**

Storage Temperature: -20 °C

TEST	Specification
Titer by Capture ELISA 2-8 ug/mL detected recombinant human TNF-alpha	Conforms
Working Dilution by Direct Immunoblotting 1-2 ug/mL detects TNF-alpha	Conforms
Working Dilution by Immunohistology 8-25 ug/mL detects human TNF-alpha in fixed human peripheral blood leukocytes	Conforms
Formulation Lyophilized from 0.2 um filtered solution of PBS (pH 7.4) with 5% Trehalose	Conforms
Bioactivity 0.02-0.08 ug/ml	Conforms

Specification: PRD.0.ZQ5.10000049026

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.