


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
trans-Styrylacetic acid - 96%

Product Number: 155322
CAS Number: 1914-58-5
MDL: MFCD00002783
Formula: C₁₀H₁₀O₂
Formula Weight: 162.19 g/mol



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Crystals
Proton NMR spectra	Conforms to Structure
Titration with NaOH	95.5 - 104.5 %
GC (area %)	≥ 95.5 %

Specification: PRD.1.ZQ5.10000050456

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.