

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Alcohol Oxidase from Hansenula sp. – vacuum-dried powder, ≥ 0.6 units/mg solid

Product Number:

A0438

CAS Number:

9073-63-6

MDL:

MFCD00130452

Storage Temperature:

-20 °C

TEST**Specification**

units/mg solid

 ≥ 0.6

One unit will oxidize 1.0 micromole of MeoH to Formaldehyde per minute at pH 7.5 at 25 deg C.

Specification: PRD.0.ZQ5.10000027317

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.