


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sulfamethazine sodium salt - ≥98%

Product Number: **S5637**
CAS Number: 1981-58-4



Formula: C₁₂H₁₃N₄NaO₂S
Formula Weight: 300.31 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
50 mg/ml H ₂ O	
Loss on Drying	≤ 6.0 %
Purity by Titration	≥ 99.0 %
Redox Titration	
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000010408

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.