


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(R)-DTBM-SEGPHOS*

Product Number: **692484**
CAS Number: 566940-03-2
MDL: MFCD09753003
Formula: C74H100O8P2
Formula Weight: 1,179.53 g/mol



TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Conforms to Requirements
Solid, Powder, Crystals, Crystalline Powder and/or Chunks	
Proton NMR spectra	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Optical Rotation (Degrees)	-55.0 - -45.0
Degrees	
C= 1% CH2Cl2	
Purity (HPLC)	≥ 94.0 %
Vendor Information	Conforms
Product of Takasago International Corporation	

Specification: PRD.1.ZQ5.10000021470

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.