

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Verrucarin A from Myrothecium sp. - ≥85.0% (HPLC)

**Product Number:** V4877 CAS Number: 3148-09-2

Formula: C27H34O9 Formula Weight: 502,55 g/mol

Storage Temperature: -20 ℃



## **TEST**

Appearance (Form)
Appearance (Colour)
Solubility (Solvent)
Solubility (Conc)
Solubility (Turbidity)
Solubility (Color)
Purity (HPLC)

Recommended Retest Period

3 Years

Specification: PRE.1.ZQ5.10000002990

## **Specification**

Pow der

White to White w/ Yellow Cast

Chloroform

9.80 - 10.20 mg/ml

Clear

Colorless to Faint Yellow

≥ 85.0 %

-----

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.