

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
LanolinProduct Number: L7387
CAS Number: 8006-54-0
MDL: MFCD00081740

TEST	Specification
Appearance (Color)	Yellow
Appearance (Form)	Semi-solid
Water (by Karl Fischer)	≤ 0.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
c = 100 mg/ml; Chloroform	
Solubility (Color)	Yellow to Yellow-Green

Specification: PRD.1.ZQ5.10000100959

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.