

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium 4-hydroxyphenyltrifluoroborate - 95%

Product Number:

738417

CAS Number:


1015082-71-9

Formula:

C₆H₅BF₃KO

Formula Weight:

200.01 g/mol

**TEST****Specification**

Appearance (Color)

White and Faint to Dark

Appearance (Form)

Yellow and Faint Brown to Dark Brown

Proton NMR spectrum

Powder or Crystals

Phosphorus NMR Spectra

Conforms to Structure

Carbon

Conforms to Structure

ICP Major Analysis

33.9 - 38.2 %

Confirms Potassium and Boron Components

Confirmed

Specification: PRD.1.ZQ5.10000008662

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.