

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-FLIP $\gamma/\delta$ , C-Terminal antibody produced in rabbit - 1 mg/mL, affinity isolated antibody, buffered aqueous solution

Product Number:

F9925

MDL:

MFCD02096351

TEST	Specification
------	---------------

Formulation Solution in 0.01M PBS and 0.02% sodium azide	Conforms
Concentration 1 mg/mL	Conforms
Working Dilution Western Blot: 1:1,000 dilution Immunocytochemistry: 2ug/mL Immunofluorescence: 10 ug/mL	Conforms

Specification: PRD.1.ZQ5.10000103081

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

