


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Chloro-2-methoxypyridine - 97%

Product Number: 642967
CAS Number: 13473-01-3
MDL: MFCD06254388
Formula: C₆H₆ClNO
Formula Weight: 143.57 g/mol



TEST	Specification
Appearance (Color)	Colorless
Infrared spectrum	Conforms to Structure
Refractive index at 20 ° C	1.524 - 1.528
Appearance (Form)	Liquid
Proton NMR spectra	Conforms to Structure
Purity (GC)	≥ 96.5 %
Titration by AgNO ₃ After Oxygen Combustion	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000017477

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.