

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Acetonitrile solution – contains 0.1 % (v/v) trifluoroacetic acid, for HPLC

Product Number:

574732



MDL:

MFCD00001878

Formula:

C₂H₃N

Formula Weight:

41.05 g/mol

TEST**Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Titration with NaOH

0.093 - 0.107 %

V% Trifluoroacetic Acid

HPLC Absorbance 254 nm

≤ 10 mAU

HPLC Absorbance 210 nm

≤ 20 mAU

HPLC SPECIFICATIONS CANNOT BE GUARANTEED FOR

VERSAFLOW (TM) PACKAGED MATERIAL.

Expiration Date Period

1 Year

Specification: PRD.3.ZQ5.10000005090

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.