

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Secoisolariciresinol diglucoside - $\geq 97\%$ (HPLC)

Product Number: S0202
CAS Number: 148244-82-0
Formula: C₃₂H₄₆O₁₆
Formula Weight: 686.70 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
1 mg/mL, H ₂ O	
Proton NMR spectrum	Conforms to Structure
Purity (HPLC)	$\geq 97\%$
Loss on Drying	$\leq 3.0\%$

Specification: PRD.0.ZQ5.10000021910

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.