

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Exendin Fragment 9-39 - ≥95% (HPLC)

**Product Number:** E7269  
CAS Number: 133514-43-9

Formula: C149H234N40O47S  
Formula Weight: 3,369 AMU

TEST	Specification
Appearance (Form)	Powder
Appearance (Color)	White
Solubility (Solvent)	Water
Solubility (Conc)	1.00 - 1.04 mg/ml
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
Peptide Content (AAA)	≥ 70.0 %
Purity (HPLC)	≥ 95.0 %
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000016375

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.