


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dimethyl sulfate - ≥99%

Product Number: 320293
CAS Number: 77-78-1
MDL: MFCD00008416
Formula: C₂H₆O₄S
Formula Weight: 126.13 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Appearance (Clarity)	Clear
Purity	≥ 99.0 %

Specification: PRD.1.ZQ5.10000028256

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.