

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Calpain 1 human - aqueous glycerol solution

Product Number: **C6108**

MDL: MFCD00217575

Storage Temperature: -70 °C

TEST	Specification
------	---------------

Purity by SDS Electrophoresis	Conforms
> /= 95%	
mg protein/ml	≥ 1
units/mg protein	≥ 1000
Enzymatic Activity	
One unit will hydrolyze 1 picomole Suc-LLVY-AMC per minute at 25 deg C.	
Tested For Infectious Agents	Tested

Specification: PRD.2.ZQ5.10000006968

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.