

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cathepsin G from human sputum - ≥ 0.2 units/mg protein, salt-free, lyophilized powder

Product Number:

C5108

CAS Number:

107200-92-0

MDL:

MFCD00130750

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

μg protein/100 μg vial (Bradford)	≥ 100
units/mg protein	≥ 0.2
One unit will hydrolyze 1.0 micromole of N-Suc-ALA-ALA-PRO-PHE-P-Nitroanilide per minute at pH 7.5 at 25 deg C.	
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000074358

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.