

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ethidium bromide solution – BioReagent, for molecular biology, 10 mg/mL in H₂O

Product Number:


E1510

MDL:

MFCD00011724

Formula Weight:

394.31 AMU



TEST

Specification

Appearance (Color)

Dark Red

Appearance (Form)

Liquid

Suitability

Suitable

Suitable for use in gel electrophoresis.

Purity (HPLC)

≥ 95 %

Recommended Retest Period

2 years

Specification: PRD.1.ZQ5.10000008597.000

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.