

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium(IV) oxide, rutile - 99.995% trace metals basis

Product Number:**204730**

CAS Number:

1317-80-2

MDL:

MFCD00011269

Formula:

O₂Ti

Formula Weight:

79.87 g/mol

**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Powder or Chunks

X-Ray Diffraction

Conforms to Structure

Loss on Ignition

≤ 1.0 %

Typically

Purity

Meets Requirements

99.995% Based on Trace Metal Analysis

Trace Metal Analysis

≤ 55.0 ppm

Specification: PRD.1.ZQ5.10000003815

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.