


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 428043
CAS Number: 9002-88-4
MDL: MFCD00084423
Formula: C₂H₄



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Pellets
Infrared spectrum	Conforms to Structure
Miscellaneous Assay	20 - 30
Melt Index: g/10 minutes	
190 Degrees Celsius, 2.16kg	
Density	0.9215 - 0.9255 g/ml
Miscellaneous Assay	320 - 480 ppm
Slip Agent	

Specification: PRD.0.ZQ5.10000030873

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.