

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
(S)-(-)-2-Methyl-CBS-oxazaborolidine solution - 1 M in THF

Product Number: **674648**

MDL: MFCD00078439

Formula: C₁₈H₂₀BNO

Formula Weight: 277.17 g/mol



TEST

Specification

Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration with NaOH	0.95 - 1.10 M
Optical Rotation	
Report Results	
C= 1%, Toluene	

Specification: PRD.1.ZQ5.10000004948

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.