

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Pomalidomide-C₉-CO₂H - ≥95%

Product Number: 901525



TEST	Specification
------	---------------

Appearance (Color)	White to Very Dark Yellow to Very Dark Brown
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
Carbon Content	59.2 - 65.4 %
Nitrogen Content	8.9 - 9.9 %
Purity (HPLC)	≥ 95.0 %
LCMS	Confirmed
444 Molecular Ion	

Specification: PRD.0.ZQ5.10000101533

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.