

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-(Dimethylamino)ethyl acrylate - contains <2,000 ppm MEHQ as inhibitor, 98%

Product Number:**330957**

CAS Number:

2439-35-2

MDL:

MFCD00038233

Formula:


C₇H₁₃NO₂

Formula Weight:

143.18 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless to Dark Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

MEHQ as Inhibitor

< 2000 ppm

Specification: PRD.0.ZQ5.10000038778

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.