

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hexylzinc bromide solution - 0.5 M in THF

Product Number:**498769**

CAS Number:

124397-96-2

MDL:

MFCD00671999

Formula:


C₆H₁₃BrZn

Formula Weight:

230.46 g/mol

Storage Temperature:

2 - 8 °C



TEST	Specification
------	---------------

Purity Titration by EDTA Complexometric	0.45 - 0.65 M
--	---------------

Specification: PRD.1.ZQ5.1000004895

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.