

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)**Product Specification**

Product Name:

4-Methyl-3-nitrobenzaldehyde 98%

Product Number:

**518328**

CAS Number:

31680-07-6

MDL:


MFCD00017011

Formula:

C<sub>8</sub>H<sub>7</sub>NO<sub>3</sub>

Formula Weight:

165.15 g/mol

**TEST****Specification**

Appearance (Color)

White to Yellow

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Carbon

56.7- 59.6 %

Nitrogen

8.2- 8.7 %

Purity (GC)

≥97.5 %

**Remarks:**

Specification Date : 06/30/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.