

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)

Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Anti-Tumor Necrosis Factor Soluble Receptor Type II antibody produced in goat – IgG fraction of antiserum,  
lyophilized powder

Product Number: T2315

MDL: MFCD01633483

TEST	Specification
Formulation	Conforms
Lyophilized from a 0.2 µm filtered solution in PBS and 5% trehalose	
Working Dilution	Conforms
By Western Blot	
1 ug/mL detects rhs TNF RII at 5 ng/lane under non-reducing conditions	
Immunocytochemistry	Conforms
Concentration of 5-15 ug/ml detects TNF RII in human PBLs.	

Specification: PRD.0.ZQ5.10000088875

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.