


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Furonitrile - 99%

Product Number: **159573**
CAS Number: 617-90-3
MDL: MFCD00003223
Formula: C₅H₃NO
Formula Weight: 93.08 g/mol

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

Colorless to Yellow to Orange to Brown

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000026038

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.