


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Bromo-2-chloropyridine - 95%

Product Number: 551902
CAS Number: 53939-30-3
MDL: MFCD01318951
Formula: C₅H₃BrClN
Formula Weight: 192.44 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow and Faint Beige to Beige
Appearance (Form) Powder, Crystals, Crystalline Powder and/or Chunks	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	94.5 - 105.5 %
Purity (GC)	≥ 94.5 %

Specification: PRD.1.ZQ5.10000033882

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.