


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-(2-Aminophenyl)indole - 97%

Product Number: 376868
CAS Number: 32566-01-1
MDL: MFCD00134432
Formula: C₁₄H₁₂N₂
Formula Weight: 208.26 g/mol



TEST

Specification

Appearance (Color)	Off-White to Yellow and Faint Beige to Beige and Faint Brown to Brown
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	96.5 - 103.5 %
Purity (TLC)	≥ 97 %
Solubility (Turbidity)	Clear
Solubility (Color)	Yellow
c = 25 mg/mL; Methanol	

Specification: PRD.0.ZQ5.10000036504

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.