


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Butyltriphenylphosphonium bromide - 99%

Product Number: B102806
CAS Number: 1779-51-7
MDL: MFCD00011855
Formula: C₂₂H₂₄PBr
Formula Weight: 399.3 g/mol



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

98.5 - 101.5 %

Specification Date : 04/29/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.