


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ponasterone A – ≥65%

Product Number: **P3490**
CAS Number: 13408-56-5
MDL: MFCD00272144
Formula: C₂₇H₄₄O₆
Formula Weight: 464.63 g/mol
Storage Temperature: -20 °C



TEST	Specification
Solubility Soluble in HOAc, EtOH, MeOH, DMSO	Conforms
Appearance (Color)	White
Appearance (Form)	Powder
EmM in EtOH	≥ 8.0
Wavelength	242 - 244 nm
Purity (HPLC)	≥ 75 %A
Purity Based on HPLC and UV combined and EmM = 12.4	≥ 65 % wt

Specification: PRD.1.ZQ5.10000045658

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.