


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Asp-Phe - 96%

Product Number: 421669
CAS Number: 13433-09-5
MDL: MFCD00063155
Formula: C₁₃H₁₆N₂O₅
Formula Weight: 280.28 g/mol



TEST

Specification

Appearance (Color)	White to Beige
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	53.2 - 58.2 %
Nitrogen	9.5 - 10.4 %
Optical Rotation C = 1%; 0.5N HCl	9.5 - 15.5 deg
Purity (HPLC)	≥ 95.5 %

Specification: PRD.0.ZQ5.10000049136

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.