

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Tumor Necrosis Factor Receptor II antibody produced in mouse - clone 22221, purified immunoglobulin, lyophilized powder

Product Number: T1815

MDL: MFCD01633483

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Lyophilized from a 0.2 µm-filtered solution in PBS with carbohydrates	
Endotoxin Assay	Conforms
< 0.1 EU/µg antibody	
Immunocytochemistry	Conforms
8-25 µg/ml recommended working dilution with immersion fixed human peripheral blood lymphocytes	

Specification: PRD.4.ZQ5.10000006465

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.