

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2-Dioleoyl-sn-glycero-3-phospho-L-serine sodium salt - ~95%

Product Number: **P1060**

MDL: MFCD00133436
Formula: C₄₂H₇₇NNaO₁₀P
Formula Weight: 810.03 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/ml, CHCl ₃ :Methanol(9:1)	
Proton NMR spectrum	Conforms to Structure
Sodium (Na)	2.1 - 3.4 %
Phosphorus (P)	3.6 - 4.0 %
Purity (TLC)	≥ 93 %
(approx 95%)	

Specification: PRD.1.ZQ5.10000011067

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.