

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

ProteoMass™ Cytochrome c MALDI-MS Standard - vial of 10 nmol, (M+H⁺) 12,361.96 Da by calculation

Product Number:

C8857

CAS Number:

9007-43-6

MDL:

MFCD00130890

TEST	Specification
------	---------------

Protein (nanomoles/vial)

8.0 - 18.0

by Micro BCA

µg/vial

98.9 - 222.5

by Micro BCA

Ave of 3 m/z measurements

12300 - 12424

for [M+ H]

Suitability

Suitable

Tested and found suitable as a standard for MALDI-MS.

Specification: PRD.1.ZQ5.10000006909

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.