


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chelerythrine chloride - ≥95% (TLC), powder

Product Number: **C2932**
CAS Number: 3895-92-9
MDL: MFCD00060717
Formula: C₂₁H₁₈ClNO₄
Formula Weight: 383.82 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	Yellow to Orange
Appearance (Form)	Powder
% Water by KF	≤ 12
Carbon (anhydrous)	62.0 - 68.7 %
Nitrogen (anhydrous)	3.3 - 4.0 %
Purity (TLC)	≥ 95 %
Solubility (Color)	Yellow to Orange
Solubility (Turbidity)	Clear
2 mg + 1 mL, DMSO	

Specification: PRD.1.ZQ5.10000043334

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.