

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)


Email USA: [techserv@sial.com](mailto:techserv@sial.com)

Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Adenosine 3',5'-cyclic monophosphate tris salt - ≥97% (HPLC), powder

**Product Number:** A3262  
**CAS Number:** 102029-77-6  
**MDL:** MFCD00058346  
**Formula:** C<sub>10</sub>H<sub>12</sub>N<sub>5</sub>O<sub>6</sub>P·C<sub>4</sub>H<sub>11</sub>NO<sub>3</sub>  
**Formula Weight:** 450.34 g/mol  
**Storage Temperature:** -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H <sub>2</sub> O	
Water (by Karl Fischer)	≤ 7.0 %
Phosphorus (P) anhydrous	6.37 - 7.17 %
EmM Anhydrous	14.6 - 16.0
At pH 7.0	
Wavelength	258 - 260 nm
Purity (HPLC)	≥ 97 %

Specification: PRD.0.ZQ5.10000029733

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.