

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Glass beads, unwashed – 212–300 µm (50–70 U.S. sieve)

Product Number:

G9143

TEST	Specification
------	---------------

Appearance (Color)

Light Grey to Blue to Green

Appearance (Form)

Beads

Mesh Size

≥ 90 %

212–300 µm

Specification: PRD.1.ZQ5.10000043845

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.