IGMA-ALDRICH

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

NiSO4 • 7H2O

Product Name: Nickel(II) sulfate heptahydrate - 99.999% trace metals basis

| Product Number: | 203890 |
|-----------------|--------------|
| CAS Number: | 10101-98-1 |
| MDL: | MFCD00149814 |
| Formula: | NiO4S · 7H2O |
| Formula Weight: | 280.86 g/mol |

| TEST | Specification |
|--|----------------------------------|
| Appearance (Color) | Blue-Green |
| Appearance (Form) | Crystals or Crystals with Chunks |
| Titration with EDTA | 20.0 - 23.0 % |
| % Ni | |
| ICP Major Analysis | Confirmed |
| Confirms Nickel and Sulfur Components | |
| Purity | Meets Requirements |
| 99.999% Based On Trace Metals Analysis | |
| Trace Metal Analysis | <pre>< 15.0 ppm</pre> |

Specification: PRD.0.ZQ5.10000062904

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.