


Product Specification

Product Name:
Amikacin sulfate salt - aminoglycoside antibiotic

Product Number: **A2324**
CAS Number: 149022-22-0



Storage Temperature: 2 - 8 °C

| TEST | Specification |
|------------------------|--------------------|
| Appearance (Color) | White to Off White |
| Appearance (Form) | Powder |
| Identification (USP) | Conforms |
| Loss on Drying | ≤ 13.0 % |
| Assay (Dry Basis) | 691 - 806 µg/mg |
| pH | 6.0 - 7.3 |
| Expiration Date Period | ----- |
| 4 Years | |

Specification: PRE.6.ZQ5.10000020763

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

