

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Vasoactive Intestinal Peptide human, porcine, rat – ≥95% (HPLC), powder

Product Number: V3628
CAS Number: 40077-57-4
MDL: MFCD00167535
Formula: C147H238N44O42S
Formula Weight: 3,325.80 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
0.1% in 1% acetic acid	
% Purity	≥ 95
% Peptide content	≥ 70

Specification: PRD.0.ZQ5.10000070835

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.