


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Amino-2-norbornanecarboxylic acid

Product Number: A7902
CAS Number: 20448-79-7
MDL: MFCD00167580
Formula: C₈H₁₃NO₂
Formula Weight: 155.2 AMU



TEST	Specification
Appearance (Color) White to Off-White	White to Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	97.0 - 103.0 %
Purity (TLC)	≥ 98 %
Solubility (Turbidity)	Clear to Hazy
Solubility (Color) Solubility Concentration: 50mg/ml in 1N NH ₄ OH (AQ)	Colorless to Light Yellow

Specification: PRD.0.ZQ5.10000016418

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.