

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(tetrahydrofuran) - average $M_n \sim 2,900$

Product Number:

345334


CAS Number:

25190-06-1

MDL:

MFCD00148879

Formula:

H(OCH₂CH₂CH₂CH₂)_nOH

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Solid
Infrared spectrum Can be run as a melt	Conforms to Structure
Molecular Weight	2826 - 3117
Hydroxyl Value (mg KOH/g)	37.0 - 40.0
Color Test	≤ 90 APHA

Specification: PRD.1.ZQ5.10000096584

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.