

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methoxyperfluorobutane - 99%, mixture of n- and iso-butyl isomers

Product Number: 464309

MDL: MFCD00674093

Formula: C5H3F9O Formula Weight: 250.07 AMU F₃C F F OCH₃ F₃C F OCH₃

TEST Specification

Appearance (Color)

Appearance (Form)
Infrared spectrum

GC Purity (sum of isomers)

Specification: PRD.0.ZQ5.10000018771

Colorless

Liquid

Conforms to Structure

> 98.5 %A