

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ethyl N-Boc-piperidine-4-carboxylate - 97%

Product Number:

665150

CAS Number:

142851-03-4

MDL:


MFCD01763998

Formula:

C₁₃H₂₃NO₄

Formula Weight:

257.33 g/mol

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Viscous Liquid

Proton NMR spectra

Conforms to Structure

Purity (TLC)

≥ 97 %

Solubility (Turbidity)

Clear to Slightly Hazy

25mg/ml, CHCl₃

Solubility (Color)

Colorless to Orange

Specification: PRD.1.ZQ5.10000033269

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.