

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ethylenediaminetetraacetic acid tripotassium salt dihydrate - Vetec™ reagent grade, 98%

Product Number:

V900077

CAS Number:


65501-24-8

Formula:

C₁₀H₁₃K₃N₂O₈·2H₂O

Formula Weight:

442.5 g/mol



TEST

Specification

Appearance (Colour)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration

97.5 - 102.5 %

(COMPLEXOMETRIC)

Solubility (Colour)

Colorless

Solubility (Turbidity)

Clear

AT 160 MG/ML OF 3 M POTASSIUM HYDROXIDE

Specification: PRE.1.ZQ5.10000002604

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.