


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methaneseleninic acid - 95%

Product Number: 541281
CAS Number: 28274-57-9
MDL: MFCD03093983
Formula: CH₄O₂Se
Formula Weight: 127.00 g/mol



TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Titration with NaOH	94.5 - 105.5 %

Specification: PRD.2.ZQ5.10000005010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.