

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Citronellyl acetate - mixture of isomers, natural, ≥95%, FG

Product Number: **W231105**

CAS Number: 150-84-5



Formula: C₁₂H₂₂O₂
Formula Weight: 198.30 g/mol

TEST

Specification

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 °C	1.440 - 1.450
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 95.0 %
Purity GC	4 - 6 %
Citronellol	
Specific Gravity	0.883 - 0.893
At 25°C	
Acid Value	≤ 1.0 ml
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
60 Months	

Specification: PRD.1.ZQ5.10000040806

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.